

Converting Raw LC-HRMS/MS Files into mzML files

B. Place

6/22/2022

In order to use the data analysis tools, all files must be converted into *.mzML files. In order to convert proprietary vendor files, all users *must* download ProteoWizard MSConvert tool, which can be acquired by going to <https://proteowizard.sourceforge.io/> and downloading the most recent version of ProteoWizard.

Once downloaded and install, follow the next steps to convert the raw file(s) to *.mzML format.

1) Run ProteoWizard MSConvert program

Select the program from the start Menu > Proteowizard > MSConvert




Figure 1: Initial screen when running MSConvert

2) Select the files to be converted

The MSConvert software can convert the following vendors:

Thermo (*.raw), Waters (*.raw), SCIEX (*.wiff2), Agilent (*.D), Shimadzu (*.LCD,) Bruker (*.D)

- Select the files using the **Browse** button.

Note: If selecting only one file at a time, you must press the *Add* button to include the file in the list.

- Select the *Output Directory* Folder

Note: It will default to the same directory as the original file.

After loading the file, the program should look like the below image.




Figure 2: MSConvert with data file loaded

3) Adjust the mzML file parameters

- Select **mzML** from the *Output format* input and make sure the *Extension* input is blank
- Select **64-bit** under *Binary encoding precision*
- Check the box next to *Write Index*, *TPP Compatibility*, and *Use zlib compression* and leave all other boxes unchecked.

See the image above for the proper selection

4) Add the conversion filters

To use the data analysis tool the following filters *must* be used, after selecting the proper parameters click the **Add** button:

Select the *Peak Picking* filter and include the following parameters:

- **Algorithm:** Vendor or CWT


Note: Vendor filter does not work for Waters files

- **MS Levels:** 1 -

If using CWT:

- **Min. SNR:** 0.1
- **Min. peak spacing:** 0.1

Parameters should look as follows



The screenshot shows a 'Filters' panel with a dropdown menu set to 'Peak Picking'. Below this, there is an 'Algorithm:' dropdown menu set to 'Vendor (does not work for UNIFI, and it MUST be the first filter!)'. Underneath, there are three input fields: 'MS Levels:' with a value of '1' and a hyphen, 'Min SNR:' with a value of '0.1', and 'Min peak spacing:' with a value of '0.1'. At the bottom of the panel are two buttons: 'Add' and 'Remove'.


Figure 3: Peak picking filter parameters

Press the **Add** Button

Select the *Threshold Peak Filter* and include the following parameters:

- **Threshold type:** Absolute intensity
- **Orientation:** Most intense
- **Value:** 1

Parameters should look as follows



The image shows a software interface window titled "Filters". Inside the window, there is a dropdown menu at the top set to "Threshold Peak Filter". Below this, there are three configuration options: "Threshold type:" with a dropdown menu set to "Absolute intensity", "Orientation:" with a dropdown menu set to "Most intense", and "Value:" with a text input field containing the number "1". At the bottom of the window, there are two buttons: "Add" and "Remove".

Figure 4: Threshold peak filter parameters

Press the **Add** Button

Note: The TitleMaker filter shows up upon opening MSConvert every time, this can be included and will not affect the data analysis.

5) Press Start Button

This will convert all added files to *.mzML format and put them in the assigned Output directory.




Figure 5: MSConvert Start button

6) Save the settings as a preset using the *Save Preset* dropdown menu (optional)




Figure 6: MSConvert presets menu